

Surface Staining with H33342 Cell Cycle Analysis on Fixed Cells

Reagents:

- Antibodies
- Fc Block
- PBS/HBSS
- PBS + 2 % FBS + 0.01% Sodium Azide
- 5 mg/mL Hoechst (H33342/33258) stock solution in PBS
- Perm Buffer: PBS + 0.5% Tween 20 + 2% FBS
- PI/7-AAD or Live/Dead Fixable stains

Protocol:

Prepare Cells:

- Collect cells and transfer approximately 1×10^6 cells/mL into a 5mL (12 x 75 mm tube)
- Wash 2x with ice cold PBS, by centrifuging at 1000 rpm, 10 mins, 4 degrees in a swinging bucket rotor centrifuge.
- Remove supernatant
 - Avoid disturbing the pellet during aspiration and take care not to aspirate to dryness
- Disrupt the pellets by flicking the tubes and resuspend in 1mL of protein free buffer (PBS/HBSS)

Live/Dead Fixable Stain**:

- Stain cells with a live/dead fixable stain compatible with chosen surface markers and H33342
 - Note that Live-Dead violet/aqua and DAPI are not compatible with H33342
 - Alternatively, PI or 7-AAD can be used with live cells at the *end* of the staining procedure.
- Add 1 μ L of reconstituted dye to the cells and incubate for 30 minutes on ice for 30 mins.
- Wash the cells with PBS and resuspend in 100 μ L of PBS

** see manufacturer's product data sheet for more details.

Surface Staining:

Block surface Fc receptors for 10 minutes at room temperature (or 15-20 minutes at 4 degrees) – DO NOT WASH

Table 1: Fc Receptor Block Reagents

Cell Type	Fc Block Source Reagents
Human	Polyclonal human IgG 10% normal human serum diluted into PBS
Mouse	Purified antibodies directed against Fc γ II/III receptor
Rat	Purified antibodies against Fc γ III receptor

Note: Most FCM Antibody reagent companies sell blocking reagents

Surface Staining with H33342 Cell Cycle Analysis on Fixed Cells

- Stain cells with directly conjugated antibodies at optimized concentrations (see antibody titration) for 15 – 20 minutes at Room temperature or 30 minutes at 4 degrees, in the dark.
 - Prepare single stain compensation controls for each fluorochrome you intend to analyze (beads or cells)
 - Reserve some unstained cells for instrument set-up.
 - Prepare fluorescent minus one gating control samples
 - Prepare matched isotype controls for each antibody-fluorophore combination (if necessary)
- Wash cells 2x with PBS and resuspend in 200 μ L of PBS + 2% FBS.

Fixation and Permeabilization:

Primary cultures with stem cell-like properties will efflux H33342 resulting in poor saturation. In this circumstance, fixation and permeabilization is necessary to achieve good DNA content resolution.

- Fix cells by adding 200 μ L of a 4% PFA solution to the cells for final concentration of 2% PFA
- Incubate on ice for 30 mins, in the dark.
- Wash cells 2X with PBS, pellet and resuspend in 500 μ L of Perm Buffer

Nuclear Staining:

- Add 2 - 5 μ L of H33342 stock solution to the appropriate tubes
 - Do not wash.
 - The amount of H33342 needed to achieve nuclear saturation may differ between cell lines and cultures – this may require titration.

Data Acquisition:

- Acquire data on an instrument equipped with a UV excitation source (355/375 nm) for H33342 excitation with detection at 450 nm.

NOTES: Other nuclear staining dyes including DAPI the FxCycles may be used with this protocol

References:

Belloc, F., Lacombe, F., Domain, P., Lopez, F., Bernard, P., Boisseau, M. R. and Reiffers, J. (1992), Intercalation of anthracyclines into living cell DNA analyzed by flow cytometry. *Cytometry*, 13: 880–885. doi:10.1002/cyto.990130811

Belloc, F., Dumain, P., Boisseau, M. R., Jalloustre, C., Reiffers, J., Bernard, P. and Lacombe, F. (1994), A flow cytometric method using Hoechst 33342 and propidium iodide for simultaneous cell cycle analysis and apoptosis determination in unfixed cells. *Cytometry*, 17: 59–65. doi:10.1002/cyto.990170108

Surface Staining with H33342 Cell Cycle Analysis on Fixed Cells

Schmid, I., Krall, W. J., Uittenbogaart, C. H., Braun, J. and Giorgi, J. V. (1992), Dead cell discrimination with 7-amino-actinomycin D in combination with dual color immunofluorescence in single laser flow cytometry. *Cytometry*, 13: 204–208. doi:10.1002/cyto.990130216

Maria Celia Jamur and Constance Oliver, (2009), Permeabilization of Cell Membranes, *Immunocytochemical Methods and Protocols*, 588 of the series *Methods in Molecular Biology* pp 63-66.