

# Direct Immunofluorescent Staining of Cell Surface Antigens

## Reagents:

- Antibodies
- FcBlock
- PBS
- Staining Buffer (PBS + 3% FBS + 0.01% Sodium Azide)

## Direct Staining in 12 x 75 mm tubes:

- Collect cells and transfer approximately  $1 \times 10^5 - 10^7$  cells/mL into a 5mL (12 x 75 mm tube)
- Wash 2 x with ice cold PBS, by centrifuging at 1000 rpm, 10 mins, 4 degrees in a centrifuge with a swinging bucket rotor.
- Remove supernatant
  - Avoid disturbing the pellet during aspiration and take care not to aspirate to dryness
- Disrupt the pellets by flicking the tubes.
- Block surface Fc receptors for 10 minutes at room temperature (or 15-20 minutes at 4 degrees) – DO NOT WASH

**Table 1:** Fc Receptor Block Reagents

Cell Type	Fc Block Source Reagents
Human	Polyclonal human IgG 10% normal human serum diluted into PBS
Mouse	Purified antibodies directed against FcγII/III receptor
Rat	Purified antibodies against FcγIII receptor

Note: Most FCM Antibody reagent companies sell blocking reagents

- Stain cells with directly conjugated antibodies at optimized concentrations (see antibody titration) for 15 – 20 minutes at RT or 30 minutes at 4 degrees, in the dark.
  - Prepare single stain compensation controls for each fluorochrome you intend to analyze (beads or cells)
  - Include an unstained sample
  - Prepare fluorescent minus one gating control samples
  - Prepare matched isotype controls for each antibody-fluorophore combination (if necessary)
- Wash cells 2x with PBS and resuspend in staining buffer.
- Add a dead cell marker (PI or 7-AAD), incubate 10-15 minutes at RT (do not wash)
  - Choice of marker depends on panel.
- Analyze cells by flow cytometry as soon as possible
  - If you cannot analyze the cells within an hour of staining, consider fixing with 2%PFA ([see protocol](#)). This will preserve the samples for up to one week.
  - Do not add 7-AAD or PI to samples that require fixation. A different choice of dead cell marker is required for fixed cells.

## Direct Staining in 96-well microtitre plates:

- Collect cells and transfer approximately  $1 \times 10^5 - 1.5 \times 10^6$  cells to each assay well.
  - Add less cells/well when working with cultured cell lines

## Direct Immunofluorescent Staining of Cell Surface Antigens

- Wash 2 x with ice cold PBS, by centrifuging at 1000 rpm, 10 mins, 4 degrees in a centrifuge with a swinging bucket rotor.
- Remove supernatant with a multi-channel pipette and disrupt the pellets by pipetting.
  - Avoid disturbing the pellet during aspiration and take care not to aspirate to dryness
- Continue as above for the remaining steps.

U of T Faculty of Medicine Flow Cytometry Facility