

Flow Cytometry of Yeast DNA

1. Remove 1 mL of cells from liquid culture at 0.5 to 1.0 OD₆₀₀ and put into an eppendorf tube (be sure to always include an asynchronous wild-type haploid culture to calibrate the flow cytometer).
2. Spin down the cells in a microfuge (14,000 rpm for 1 minute).
3. Aspirate off the supernatant (be careful not to suck up pellet!).
4. Add 1 mL of 70% EtOH and resuspend the pellet by vortexing.
5. Let the cells sit at room temperature for at least 1 hour, then store at 4°C or continue with the protocol.
6. Spin down 0.5 mL of cells and carefully aspirate off the supernatant. Keep the other 0.5 mL of ethanol fixed cells at 4°C.
7. Resuspend the cells in 0.5 mL ddH₂O.
8. Spin down the cells and aspirate off the supernatant (the pellet may not be visible, so be careful when aspirating).
9. Resuspend in 200 µL of RNase A solution (4 µL of 10mg/mL RNase A + 196 µL 50 mM Tris-Cl pH 8.0; make up a mix for all your samples).
10. Incubate the cells at 37°C for 2-4 hours.
11. Spin down the cells and aspirate off the supernatant.
12. Resuspend the pellet in 200 µL proteinase K solution (2 mg/mL proteinase K in 50 mM Tris-Cl pH7.5; make up a mix for all your samples).
13. Incubate at 50°C for 30-60 minutes.
14. Spin down the cells and aspirate off the supernatant.
15. Resuspend in 200-400 µL FACS buffer (can leave in FACS buffer at 4°C, but no longer than 1 week).
16. Transfer 10 µL of cells into 96-well plate. Add 200 µL SybrGreen (diluted 5,000X from stock in 50 mM Tris-Cl pH7.5) to each well. Samples are light sensitive, keep in dark place where possible.
27. Sonicate each sample for 3 seconds at 10% power (5-10 watt output).
28. Run samples immediately on the Flow Cytometer.

Materials

70% Ethanol (250 mL):

Mix 184.2 mL of 95% ethanol with 65.8 mL H₂O

Filter sterilize

1M Tris-Cl pH 7.5/8.0 (500 mL each):

Add 60.57g of Tris into 300 mL H₂O

Calibrate pH meter

Add HCl until solution reaches appropriate pH

Bring volume to 500 mL with H₂O

Autoclave

5M NaCl (500 mL):

Add 146.1g of NaCl to 300 mL H₂O

Bring solution to near 500 mL and add heat to dissolve NaCl

Bring volume to 500 mL with H₂O

Autoclave

1M MgCl₂ (100mL):

Add 20.33g MgCl₂ to H₂O

When dissolved, bring to 100 mL with H₂O

Autoclave

RNAse A Solution:

10 mg/mL of RNase A stored at -20°C freezer

Working Solutions:

50mM Tris-Cl pH 7.5/8.0 (250mL):

Add 12.5 mL of 1M stock to 237.5 mL of H₂O

Filter sterilize

FACS Buffer (100mL):

200 mM Tris-Cl pH 7.5	20 mL of 1M stock
200 mM NaCl	4 mL of 5M stock
78 mM MgCl ₂	7.8 mL of 1M stock
H ₂ O	68.2 mL

Reference:

Dunham, M.J., Gartenberg, M.R., and Brown, G.W. (2015). *Methods in Yeast Genetics and Genomics, A Cold Spring Harbor Laboratory Course Manual, 2015 Edition*. Cold Spring Harbor, NY: Cold Spring Harbor Laboratory Press.