

# BD FACSymphony™ A3-II



RED 637	780/60	730/45	670/30					
	APC-Cy7, Zombie NIR, Live/Dead NIR	Alexa 700, NovaRed710	APC, AF647, NovaRed660					
YELLOW-GREEN 561	780/60	710/50	670/30	610/20	586/15			
	PE-Cy7	PE-Cy5.5, NovaYellow700	PE-Cy5, 7-AAD, NovaYellow660	PE-CF594, PE-e610, mCherry, PI, NovaYellow610, PE-Dazzle	PE, tdTomato, RY586			
BLUE 488	810/40	710/50	670/30	610/20	530/30			
	BB790, BB750, PB 747	BB700, PerCP-e710, PB 704	BB660, PerCP-Cy5.5, NovaBlue-120S	BB630, NovaBlue610	BB515, CFSE, FITC, Alexa 488, GFP, NovaBlue530			
VIOLET 405	810/40	750/30	710/50	677/20	605/40	586/15	525/50	431/28
	BV786	BV750	BV711	BV650	BV605, SB600	Live/Dead Yellow, BV570	BV510, BV480, CellTracker Violet, Live/Dead aqua, Zombie Aqua, VioGreen, CFP	BV421, Pacific Blue, Zombie Violet, eFluoro450, TagBFP, Alexa 405, VioBlue
ULTRAVIOLET 349	810/40	735/30	670/25	610/20	580/20	515/30	379/28	
	BUV805	BUV737	BUV661	BUV615	BUV563	BUV496, Hoechst, Zombie UV, DAPI, FixableBlue	BUV395	

# BD FACSymphony™ A3-II: Relative Fluorochrome Ranking Measure by CD4 Stain Index

Staining Index	Laser				
High	Ultraviolet	Violet	Blue	Yellow-Green	Red
				PECy7	
				PE-CF594	
				PECy5	
		BV421	BB630-P2, BB515		
			BB700	PE	
			BB790, BB66-P2	RY586	APC-R700, R718
		BV750, BV786			
	BUV661, BUV737	BV711, BV605			
	BUV805				APC
	BUV615	BV650			APC-CY7
	BUV395	BV480			AF700
		BV510			
		BV570			
Low	BUV563, BUV496		PerCP-Cy5.5, FITC		

